Project: Sapien 190 Date: 10.05.2010 Version: 0.01

Product architecture document for Sapien 190

Embedded Real-Time Systems (TI-IRTS) Spring 2010

Peter Høgh Mikkelsen (20087291) Anders Block Arnfast (20085515) Kim Bjerge (20097553)

Document history

Date	Version	Description	Author
24.04.2010	0.00	Initial Version	KBE
10.05.2010	0.01	Updated with Use Case View and initial UML diagrams for first delivery UC#1	KBE

Template User guide:

- To be deleted in the final document

This Template is a highly modified version of a RUP template (Rational Unified Process) for a software architecture document. The modifications includes elements from the Danish Structured Development Handbooks (SPU) SW documentation guide. The template is based on the famous "4+1" View Model described by Philippe Kruchten, "The 4+1 View Model of Architecture", IEEE Software Nov.1995.

The Template is written as a **Microsoft Word** .dot file. A double click with the mouse on the xx.dot file in the file manager will open a new word document with the default name: Document 1, including the information in the template.

Start with modifying the document title and subject. Document title and subject (e.g. System/product <name>) is modified in **File/Properties** – **Summary**. Title and subject will automatically be inserted on the front page and in the page header. Save the document with a new name.

If the .dot file is installed in the directory where word looks for template files, it will be possible to activate a new document based on this template by selecting **File/New** from the Word program.

For helping the user with writing the document, help text is included for each section in the document. This help text is written as "Hidden Text" in Word.

To show (enable) the help text on the screen:

Menu choice: Tools/Options - View

and check the check box: Hidden text (alternatively remove the mark – for removing the help text).

To print the template file with help text:

Menu choice: Tools/Options - Print

and check the check box Hidden text.

To modify the template:

Start Word and open the template by selecting: File/Open

Version 0.92 27.11.2002 side 0 af 13

Table of Contents

1. INTRODUCTION	3
1.1 Purpose and Scope	
1.3 Definitions and acronyms	
1.4 Document structure and reading guide	
1.5 Document role in an iterative development process	
2. SYSTEM OVEVIEW	
2.1 System context	
2.1 System context 2.2 System introduction	
3. SYSTEMET INTERFACES	5
3.1 Interface to human actors	5
3.2 Interface to external system actors	
3.3 Interface to hardware actors	
3.3.1 Pulse interface	5
3.3.2 Infusion pump interface	
3.4 Interface to external software actors	5
4. USE CASE VIEW	5
4.1 Overview of architecture significant Use cases	6
4.2 Use case #1: Execute and Control Simulation scenarios	
4.2.1 Use case goal	6
4.2.2 Use case scenarios	
4.3 Use case # 2: Select and Initiate Scenario scenarios	7
5. LOGICAL VIEW	8
5.1 Overview	8
5.2 Architecturally significant design packages	9
5.2.1 Continuous Package	
5.2.2 Discrete	
5.2.3 Com	
5.2.4 AbstractHW	
5.2.5 AbstractOS	
5.2.6 Application Helper Classes	
5.3.1 Use case #1: Execute and Control Simulation scenarios	
5.3.2 Use case 2. realization	
6. PROCESS/TASK VIEW	
6.1 Process/task overview	
6.3 Process/task communication and synchronization	
6.4 Process group 1.	
6.4.1 Process communication in group 1	
6.4.2 Process 1. description	
6.4.3 Process 2. description	
6.5 Process group 2.	22

7. DEPLOYMENT VIEW	22
7.1 System configurations overview	22
7.2 System configurations	
7.2.1 Configuration 1	22
7.2.2 Configuration 2.	
7.3 Node descriptions	
7.3.1 Node 1. description	
7.3.2 Node 2. description	22
8. IMPLEMENTATION VIEW	23
8.1 Overview	23
8.2 Component descriptions	23
8.2.1 Component 1	23
8.2.2 Component 2	23
9. DATA VIEW	23
9.1 Data model	23
9.2 Implementation of persistence	
10. GENEREL DESIGN DECISIONS	
10.1 Architectural goals and constraints	23
10.2 Architectural patterns	
10.3 General user interface design rules	
10.4 Exception and error handling	
10.5 Implementation languages and tools	
10.6 Implementation libraries	
11. SIZE AND PERFORMANCE	24
12. QUALITY	24
13. COMPILATION AND LINKING	24
13.1 Compilation-hardware	24
13.2 Compilation-software	
13.3 Compilation and linking process	24
14. INSTALLATION AND EXECUTING	24
14.1 Installation	24
14.2 Executing-hardware	
14.3 Executing-software	24
14.4 Execution-control (start, stop and restart)	25
14.5 Error messages	25
15 APPENDICES	25

1. INTRODUCTION

1.1 Purpose and Scope

The Sapien 190 is a patient simulator, simulating human physiological behaviour according to different patient scenarios. Scenarios are written in an open format and can be downloaded to the Sapien 190 that contains patient records from the PhysioBand database.

In conjunction with the Sapien 110 human doll, the Sapien system provides a complete simulated human interface, with EKG, ECG and respiratory measuring spots and medicine injection spots.

1.2 References

- [1] Erich Gamma et al., Design Patterns: Elements of Reusable Object-Oriented Software, Addison Wesley (GoF)
- [2] Bruce Powell Douglass, Real-Time Design Patterns: Robust Scalable Architecture for Real-Time Systems
- [3] PhysioNet and PhysioBank the research resource for complex physiologic signals http://www.physionet.org/
- [4] Project specification for PSIMU: Patient Simulator System http://kurser.iha.dk/eit/tiirts/Projekter/PSIMU-project.doc
- [5] Project specification for LMON: Local Monitor System http://kurser.iha.dk/eit/tiirts/Projekter/LMON-project.doc
- [6] Project specification for IPUMP: Infusion Pump System http://kurser.iha.dk/eit/tiirts/Projekter/IPUMP-project.doc
- [7] Project Interface Specification
 http://kurser.iha.dk/eit/tiirts/Projekter/ProjectInterfaces.doc
- [8] Requirement specification for Sapien 190

1.3 Definitions and acronyms

- **Model** A model that represents one physical individual. The model may consist of sub-models for different body subsystems. The model uses an algorithm to compute the output signals based on the input signals or patient records.
- **Model Parameters** Parameters used in the model.
- **Record** A patient record file taken from the PhysioBank database
- **Scenario** Model parameters and a collection of records taken from the PhysioBank database.
- **Scenario Configuration** A set of files that represents model parameters and patient records.
- **Signals** Signals in form of waveform files or input from external equipment
- **Simulation** A continuous mode, where the physiological output signals are updated according to the model and the scenario applied to it.
- **ECG** Electrocardiogram
- EDR ECG–Derived Respiration

Version 0.92 27.11.2002 side 3 af 13

1.4 Document structure and reading guide

Chapter 2 and 3 gives an overview of the product and the interfaces to the patient simulator in terms of other devices and user operation.

The document describes the design using the "4+1" view. For each iteration as specified in ROPES [2] we have selected one or more use cases that is used in describing the Use Case, Logical, Process, Deployment and Implementations Views. These views are described in chapters 5-8.


1.5 Document role in an iterative development process

This document is updated in using the ROPES development process [2]. The document is updated for every design cycle of the ROPES spiral microcycle. Every design microcycle is using the 4+1 view

2. SYSTEM OVEVIEW

The Sapien 190 is a compact unit with a graphical user interface (GUI) and interfaces to a range of body monitoring and injection equipment. The Sapien 190 can be operated by the instructor and connected to a bedside monitor that used the output signals from the patient simulator. Optional a medicine infusion pump can be connected to the simulator.

2.1 System context



2.2 System introduction

The instructor will be able to control the simulation like forcing a heart attack and monitor the waveform of the signals that is send to the bedside monitor. Two analogue outputs signals can be connected to a bedside monitor to display ECG and EDR signals. It is possible remotely using an Ethernet connection to the patient simulator to update and delete files in the scenarios database on Sapien 190. The optional infusion pump injects medicine into the patient and the flow of medicine can be monitored on the Sapien 190 LCD display.

Version 0.92 27.11.2002 side 4 af 13

3. SYSTEMET INTERFACES

3.1 Interface to human actors

Instructor controls the patient simulation by inputs to the LCD touch screen....

3.2 Interface to external system actors

Remotely it is possible to update the scenarios database by use of a ftp connection with a Ethernet connected to the patient simulator.

3.3 Interface to hardware actors

Analogue outputs signals to the LMON has a resolutions of 12 bits sampled at 250 hz. The analogue output has a voltage range of 0 - 4.096 volts.

3.3.1 Pulse interface

The pulse is send to a connected monitor using a RS232 connection (115200 baud, 8 databits, even parity, 1 stop bit).

Format of the pulse signal is a maximum 3 bytes ASCII value (beats/minute). It is transmitted from the PSIMU to the LMON with a rate of 1 second. The pulse value is terminated with a <CR>. <Pulse> <CR>

Example a terminal can be connected to display the pulse values updated each second: 103

107

3.3.2 Infusion pump interface

The infusion pump is connected using a RS232 connection with the protocol described below.

A fixed 64 bytes ASCII PDU is transmitted from IPUMP to PSIMU.

This PDU is transmitted every second when the pump is started.

PDU format:

4 bytes startframe (##?*)

46 bytes medicin name

12 bytes volume infused (since started)

2 bytes CRC checksum

3.4 Interface to external software actors

4. USE CASE VIEW

In the first iteration the UC #1 has been selected since it provides the main functionality and is the essential use case for the architecture of the patient simulator.

Execute and Control Simulation

Version 0.92 27.11.2002 side 5 af 13

When the patient simulation is running it will perform reading of the digitized physiologic signals and send theses "real time" values as analogue signals to local connected bedside monitoring equipments. Up to 2 analogue channels with different signals is possible to be simulated simultaneously. The simulated signals can be ECG or EDR. The pulse will be signaled to the bedside monitoring equipment as a digital signal.

4.1 Overview of architecture significant Use cases

The UC #1 has been selected for the first iteration of the ROPES spiral microcycle [2] in making the architectural, mechanistic and detailed design. This use case is significant and provides the central functionality of the patient simulator. This use case provides the basis functionality that allows the monitor to be connected being able to display the ECG and EDR signals. It also reduces the risk for developing the patient monitor since it covers all the unknown technologies of the product like:

- Reading the patient record files on the target (Linux, WFDB and target)
- Generating the analogue output signals (Writing to drivers)
- Display of signal waveform using Qt on target (Working with Qt on target)

4.2 Use case #1: Execute and Control Simulation scenarios

4.2.1 Use case goal

To simulate signals from the patient based on the selected scenario.

The waveform of signals to monitor and the status of patient are displayed.

The instructor must be able to monitor the simulated patient.

4.2.2 Use case scenarios

Scenario 1 - normal:

- 1. Opens scenario file
- 2. Search for record file
- 3. Opens record file
- 4. Continues to reads samples from record file and performs:
 - a. Use patient model to generate ECG signal
 - b. Use patient model and ECG signal to calculate EDR signal
 - c. Use patient model and ECG signal to calculate pulse
 - d. Update output signals: Pulse, ECG and EDR

Scenario 2 – normal with alternative scenario record:

- 1. Opens scenario file
- 2. Search for record file
- 3. Opens record file
- 4. Continues to reads samples from record file and performs:
 - a. Use patient model to generate ECG signal
 - b. Use patient model and ECG signal to calculate EDR signal
 - c. Use patient model and ECG signal to calculate pulse

Version 0.92 27.11.2002 side 6 af 13

- d. Update output signals: Pulse, ECG and EDR
- 5. Open alternative record file after specified simulation time and continues at 2.

Scenario 3 – normal with IPUMP:

- 1. Opens scenario file
- 2. Search for record file
- 3. Opens record file
- 4. Continues to reads samples from record file and performs:
 - a. Reads medicine and volume from IPUMP
 - b. Use patient model to generate ECG signal
 - c. Use patient model and ECG signal to calculate EDR signal
 - d. Use patient model and ECG signal to calculate pulse
 - e. Update output signals: Pulse, ECG and EDR

Scenario 4 – error opening record file:

- 1. Opens scenario file
- 2. Failed to search and open record file
- 3. Error message on LCD display

4.3 Use case # 2: Select and Initiate Scenario scenarios


Initiate the simulation using a scenario configuration chosen by the instructor.

This use case is not part of the first iteration and will be updated in the next iteration.

Version 0.92 27.11.2002 side 7 af 13

5. LOGICAL VIEW

5.1 Overview



Version 0.92 27.11.2002 side 8 af 13

5.2 Architecturally significant design packages

5.2.1 Continuous Package




Figure 1 Major Classes in Continuous Package

Version 0.92 27.11.2002 side 9 af 13




Figure 2 Boundary Classes for RealTimeThread




Figure 3 Proxy and Iterator Pattern used to access Records from the PatientModel

Version 0.92 27.11.2002 side 10 af 13




Figure 4 Record and Iterator




Figure 5 Strategy, Filter and Pipes Pattern used for PhysioModel

Version 0.92 27.11.2002 side 11 af 13




Figure 6 Strategy for PatientModel




Figure 7 Filter Pattern for PhysioModel

Version 0.92 27.11.2002 side 12 af 13




Figure 8 Façade Pattern used for interface to the Continuous Package

5.2.2 Discrete




Figure 9 Application model for discrete package

Version 0.92 27.11.2002 side 13 af 13

5.2.3 Com




Figure 10 Communication package Serial protocol

5.2.4 AbstractHW




Figure 11 Hardware Abstraction

Version 0.92 27.11.2002 side 14 af 13

5.2.5 AbstractOS




Figure 12 Abstract OS (Linux)

5.2.6 Application Helper Classes




Figure 13 Application Helper classes

Version 0.92 27.11.2002 side 15 af 13

5.3 Use case realizations

5.3.1 Use case #1: Execute and Control Simulation scenarios




Figure 14 Sequence Diagram: Create Realtime Simulator

Version 0.92 27.11.2002 side 16 af 13




Figure 15 Sequence Diagram: Create Simulation Record and Set Medicine

Version 0.92 27.11.2002 side 17 af 13




Figure 16 Sequence Diagram: Start Realtime Simulation

Version 0.92 27.11.2002 side 18 af 13

5.3.1.1 Test setup for #UC 1 continuous package SimulatorRealtime




Figure 17 Test setup for UC#1 – test setup in Rhapsody only

Version 0.92 27.11.2002 side 19 af 13





Figure 18 State machine test of UC#1 including simulation of Threads in Rhapsody

5.3.2 Use case 2. realization

Version 0.92 27.11.2002 side 20 af 13

6. PROCESS/TASK VIEW

6.1 Process/task overview



6.2 Process/task implementation

6.3 Process/task communication and synchronization

6.4 Process group 1.

6.4.1 Process communication in group 1

6.4.2 Process 1. description


6.4.3 Process 2. description

Version 0.92 27.11.2002 side 21 af 13

6.5 Process group 2.

7. DEPLOYMENT VIEW

7.1 System configurations overview



7.2 System configurations

7.2.1 Configuration 1.

7.2.2 Configuration 2.

7.3 Node descriptions

7.3.1 Node 1. description

7.3.2 Node 2. description

Version 0.92 27.11.2002 side 22 af 13

8. IMPLEMENTATION VIEW

- 8.1 Overview
- 8.2 Component descriptions
- **8.2.1** Component 1
- **8.2.2** Component 2
- 9. DATA VIEW
- 9.1 Data model
- 9.2 Implementation of persistence

10. GENEREL DESIGN DECISIONS

10.1 Architectural goals and constraints

The architectural goals is to make the system modifiable and provide high performance. Since the system is a real time system, there are constraints on how modular system can be, since modular and performance work in the opposite direction. It should be easy to extend the system new types of medicine, different kind of inputs and output. The system should be encapsulated so it's easy to port the system to different platform.

10.2 Architectural patterns

Observer Pattern:

Observer pattern help out with mass distribution of information in a system. If a lot process are interested in the same information, observer pattern is the pattern to choose. In Sapien observer pattern is used to distribute information generated by our simulation to the user interface. This a very common way to use observer pattern since it also notify the user interface when new data has arrived so it can be refreshed.

Five Layered Architecture

Layered pattern organizes domains into a hierarchical organisation based their level of abstraction. The advance of layered architectural is that it make it easier to find relevant code, and make it easier to replace whole layers, for instance if you want to port the system to another device. There used open layered architectural so it's allowed to call more than one layer down.

Version 0.92 27.11.2002 side 23 af 13

10.3 General user interface design rules

10.4 Exception and error handling

10.5 Implementation languages and tools

Implementation languas:

- C++

Tools

- Ecliplse
- TrollTech QT Creator

10.6 Implementation libraries

QT-Everywhere-4.6.1

wfdb

11. SIZE AND PERFORMANCE

12. QUALITY

13. COMPILATION AND LINKING

- 13.1 Compilation-hardware
- 13.2 Compilation-software
- 13.3 Compilation and linking process

14. INSTALLATION AND EXECUTING

- 14.1 Installation
- 14.2 Executing-hardware

14.3 Executing-software

Version 0.92 27.11.2002 side 24 af 13

- 14.4 Execution-control (start, stop and restart)
- 14.5 Error messages
- 15. APPENDICES

Version 0.92 27.11.2002 side 25 af 13